










Renú Transoral Needle Symbols Glossary

Symbol	Symbol Title	Symbol Reference #	Symbol Meaning/ Explanatory Text	Standard Number	Standard Title
	Manufacturer	5.1.1	Indicates the medical device manufacturer, as defined in EU directives 90/385/EEC, 93/42/EEC and 98/79/EC.	EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements
	Authorized Representative in the European Community	5.1.2	Indicates the Authorized Representative in the European Community.	EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements
	Use-by Date	5.1.4	Indicates the date after which the medical device is not to be used.	EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements
	Batch Code	5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified.	EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements
	Catalogue Number	5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified.	EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements
	Sterilized using irradiation	5.2.4	Indicates a medical device that has been sterilized using irradiation.	EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements
	Do not re-sterilize	5.2.6	Indicates a medical device that is not to be re-sterilized.	EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements
	Do not use if package is damaged	5.2.8	Indicates a medical device that should not be used if the package has been damaged or opened.	EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements
	Do not re-use	5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements