












Renú Transoral Needle Symbols Glossary

| Symbol | Symbol Title | Symbol Reference # | Symbol Meaning/Explanatory Text | Standard Number | Standard Title |
|---|---|--------------------|--|---------------------|---|
|  | Manufacturer | 5.1.1 | Indicates the medical device manufacturer, as defined in EU directives 90/385/EEC, 93/42/EEC and 98/79/EC. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Authorized Representative in the European Community | 5.1.2 | Indicates the Authorized Representative in the European Community. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Use-by Date | 5.1.4 | Indicates the date after which the medical device is not to be used. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Batch Code | 5.1.5 | Indicates the manufacturer's batch code so that the batch or lot can be identified. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Catalogue Number | 5.1.6 | Indicates the manufacturer's catalogue number so that the medical device can be identified. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Sterilized using irradiation | 5.2.4 | Indicates a medical device that has been sterilized using irradiation. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Do not re-sterilize | 5.2.6 | Indicates a medical device that is not to be re-sterilized. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Do not use if package is damaged | 5.2.8 | Indicates a medical device that should not be used if the package has been damaged or opened. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Do not re-use | 5.4.2 | Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Consult instructions for use | 5.4.3 | Indicates the need for the user to consult the instruction for use. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |
|  | Caution | 5.4.4 | Indicates the need for the user to consult the instruction for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. | EN ISO 15223-1:2016 | Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements |